



HFMED®

Fully automatic alternating air bed



Heat sealing technology does not leak

Removable air strips, restarting pipelines are made by heat sealing machines from various countries, with strong sealing performance, each airbag can be disassembled and cleaned separately

### WHY CHOOSE US



#### OUR ADVANTAGE

The company keeps to the concept of being oriented by science and technology and taking quality as life. Now the company has passed the certification of ISO9001:2008, ISO13485:2003, CE and SGS.



#### MANAGEMENT

"Quality first, Credit first, Service first" is our tent. Hui Feng are looking forward to sincere cooperation with you. Let's advance and pursue a wonderful future by mutual efforts.



#### SALES MARKET

all has its domestic sales amount coming out at the top and some products have ready sale in Southeast Asia. Its products can be found in all the domestic markets in China and also are sold abroad.



Official QR code

HFMED®

YOUR IDEAL SUPPLICAL LIGHT SUPPLIER

SHANGHAI HUIFENG MEDICAL INSTRUMENT CO.,LTD

Address: NO 885, Qluxing Road, Nicheng Town, Pudong New Area, shanghai, China.

Tel: +862158084992

Mob: +862158085105

Email: lisa@hui Feng-co.com

Web: www.sh-huifeng.com

HFMED®



CE

Fully automatic alternating air bed



HFMED

SHANGHAI HUIFENG MEDICAL INSTRUMENT CO.,LTD

## Product Description

### Alternate Type A Microporous Breathable

Doctors advise bedridden patients to turn over manually every 2 hours. Air mattresses cycle every 10 minutes, equal to 6 turns per hour.



Easy to replace


### 24 independent air bags are removable

Easy to clean, reduce bacteria and facilitate maintenance and replacement. The alternating air cushion uses the air ejected to take away the moisture outside the human body, thereby drying the skin, and the principle of dispersing the patient's body pressure on the mattress to effectively prevent the occurrence of bedsores.

## Main Technical Data

### Fully automatic alternating air bed

Name	Air mattress alternate type A
Material	Nylon TPU
Size	197x65x14cm
Weight	1.12kg
Uses	hospital beds, nursing home beds, etc.



### 24 gas upgrades Circulating wave jet, breathable and comfortable

Using microcomputer control to adjust the speed of inflation and deflation, the knob of the air pump is adjusted to the 123-position fluctuating state, and the air pump is double-mouthed.

Fully automatic  
alternating air bed